

Philips Sonicare AirFloss
Ultra
Interdental cleaner

w/ 1 nozzle & 1 brush head
w/ Plaque control toothbrush

AirFloss Ultra



HX8392/43

Healthier gums in 2 weeks, guaranteed*

Designed for inconsistent flossers

For those who don't floss consistently, AirFloss Ultra is the easiest way to effectively clean between teeth. AirFloss Ultra can be used with mouthwash or water and is clinically proven as effective as floss for gum health.**

Improves oral health

- Up to 99.9% plaque removal***
- Clinically proven as effective as floss for gum health**
- Helps prevent cavities between teeth

Superior technology

- Air and micro-droplet technology
- High performance nozzle
- Triple burst customizable settings

The easiest way to effectively clean between teeth

- An easy way to start a healthy habit
- Takes just 60 seconds: Point, Press, Clean!
- Fresh and clean with mouthwash

PHILIPS
sonicare


Highlights

Up to 99.9% plaque removal




Removes up to 99.9% of plaque from treated areas.**

Micro-droplet technology




Our clinically proven results are possible from our unique technology that combines air and mouth rinse or water to powerfully yet gently clean between teeth and along the gumline.

Easily start a healthy habit



Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit.

Better gum health



Philips Sonicare AirFloss Ultra is clinically proven to improve gum health as much as floss.** Helps improve gum health in as little as two weeks.

Cleans in just 60 seconds



With the AirFloss Ultra, cleaning your entire mouth takes less than 60 seconds a day. Simply select your burst frequency (single, double or triple), and hold down the activation button for continuous automatic bursts or press and release for manual burst mode.

Fresh and clean with mouthwash



Use with mouthwash for the ultimate fresh experience and anti-microbial benefits.

Helps prevent cavities




By gently bursting away plaque that brushing missed, Philips Sonicare AirFloss Ultra helps prevent cavities from forming in the spaces between your teeth.

High performance nozzle



The new Airfloss Ultra high performance nozzle amplifies the power of our air and micro-droplet technology to be more effective and efficient than ever.

Triple burst customization



AirFloss Ultra's powerful spray can be customized to your liking, with single, double or triple bursts for each press of the activation button.

Specifications

Items included

- AirFloss Ultra handle: 1
- Nozzle: 1
- Rechargeable toothbrush handle: 2 Series plaque control
- Brush heads: 1 plaque control brush head
- Charger: 1

Design and finishing

- Color: Gray

Cleaning performance

- Clinically proven results: As effective as floss, for gum health improvement, for reducing gum bleeding, for plaque reduction
- Performance: Removes up to 99.9% of plaque
- Health benefits: Cleans entire mouth in 60 sec, Improves gum health in 2 weeks
- Customizable settings: Auto-burst or manual

- burst, Single, double, triple burst
- For best results: Change nozzle every 6 months

Ease of use

- Handle: Slim ergonomic design
- Battery indicator: Handle LEDs show charge level
- Battery Life: Up to 2 weeks
- Nozzle system: Easily snaps on and off

Technical specifications

- Battery: Rechargeable
- Battery type: Lithium ION
- Charge time: 24 hours to full charge

Power

- Voltage: Multi-voltage charger

Service

- Warranty: 2-year limited warranty



Issue date 2017-06-02

Version: 6.1.1

EAN: 08 71010 37151 77

© 2017 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* or your money back

** when used in conjunction with a manual toothbrush and anti-microbial mouth rinse in patients with mild to moderate gingivitis; AirFloss is designed to help inconsistent flossers develop a healthy daily interdental cleaning routine. Please see Q&A under Support tab for further details.

*** In a lab study, actual in-mouth results may vary.