




Philips Sonicare AirFloss
Ultra
Rechargeable powered
interdental cleaner

3 nozzles

AirFloss Ultra



HX8433/05

Healthier gums in 2 weeks, guaranteed*

Designed for inconsistent flossers

For those who don't floss consistently, AirFloss Ultra is the easiest way to effectively clean between teeth. AirFloss Ultra can be used with mouthwash or water and is clinically proven as effective as floss for gum health.**

Improves oral health

- Up to 99.9% plaque removal***
- Clinically proven as effective as floss for gum health**
- Helps prevent cavities between teeth

Superior technology

- Air and micro-droplet technology
- High performance nozzle
- Triple burst customizable settings


The easiest way to effectively clean between teeth

- Total daily confidence of a fresh and clean mouth
- An easy way to start a healthy habit
- Takes just 60 seconds: Point, Press, Clean!

PHILIPS
sonicare

Highlights

Better gum health



Philips Sonicare AirFloss Ultra is clinically proven to improve gum health as much as floss.** Helps improve gum health in as little as two weeks.

Up to 99.9% plaque removal




Removes up to 99.9% of plaque from treated areas.***

Helps prevent cavities



By gently bursting away plaque that brushing missed, Philips Sonicare AirFloss Ultra helps prevent cavities from forming in the spaces between your teeth.

Micro-droplet technology



Our clinically-proven results are possible from our unique technology that combines air and mouthwash or water to powerfully yet gently clean between teeth and along the gumline.

High performance nozzle



The new Airfloss Ultra high performance nozzle amplifies the power of our air and micro-droplet technology to be more effective and efficient than ever.

Triple burst customization




AirFloss Ultra's powerful spray can be customized to your liking, with single, double or triple bursts for each press of the activation button.

Cleans in just 60 seconds




With the AirFloss Ultra, cleaning your entire mouth takes less than 60 seconds a day. Simply select your burst frequency (single, double or triple), and hold down the activation button for continuous automatic bursts or press and release for manual burst mode.

Fresh & clean with mouthwash



Fill the reservoir on the handle with either mouthwash or water, then point and press. Use with mouthwash for the ultimate fresh experience and anti-microbial benefits.

Easily start a healthy habit



Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit.

Specifications

Items included

- AirFloss Ultra handle: 1
- Nozzle: 3
- Nozzle holder: 1
- Charger: 1

Cleaning performance

- Clinically proven results: As effective as floss, for gum health improvement, for reducing gum bleeding, for plaque reduction
- Performance: Removes up to 99.9% of plaque
- Health benefits: Cleans entire mouth in 60 sec, Improves gum health in 2 weeks
- Customizable settings: Auto-burst or manual burst, Single, double, triple burst
- For best results: Change nozzle every 6 months

Ease of use

- Handle: Slim ergonomic design

- Battery indicator: Handle LEDs show charge level
- Battery Life: Up to 2 weeks
- Nozzle system: Easily snaps on and off

Technical specifications

- Battery: Rechargeable
- Battery type: Lithium ION
- Charge time: 24 hours to full charge

Power

- Voltage: Multi-voltage charger

Service

- Warranty: 2-year limited warranty

Design and finishing

- Color: Black



Issue date 2017-08-10

Version: 2.0.1

© 2017 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* or your money back

** when used in conjunction with a manual toothbrush and anti-microbial mouth wash in patients with mild to moderate gingivitis; AirFloss is designed to help inconsistent flossers develop a healthy daily interdental cleaning routine. Please see Q&A under Support tab for further details.

*** In a lab study, actual in-mouth results may vary.