

Register your product and get support at
www.philips.com/welcome

HX9383, HX9343, HX9333

Sonicare DiamondClean

Quick start guide

Kurzanleitung

Guide de mise en route

Guida di avvio rapido

Snelstartgids



Please read chapter 'Important' in the user manual before you use the appliance.

Lesen Sie bitte in der Bedienungsanleitung das Kapitel "Wichtig", bevor Sie das Gerät verwenden.


Veuillez lire le chapitre « Important » du mode d'emploi avant d'utiliser l'appareil.

Leggete il capitolo 'Importante' nel manuale per l'utente prima di iniziare a utilizzare l'apparecchio.

Lees hoofdstuk 'Belangrijk' in de gebruiksaanwijzing voordat u het apparaat gebruikt.

PHILIPS

CONTENT / INHALT / SOMMAIRE / SOMMARIO / INHOUD



A
Brush head
Bürstenkopf
Tête de brosse
Testina
Opzetborstel

B
Handle
Handstück
Manche
Impugnatura
Handvat


C
Charger glass
Glas des Ladegeräts
Verre de charge
Bicchieri per ricarica
Oplaadglas

D
Charging stand
Ladegerät
Base de recharge
Supporto di ricarica
Oplaadstandaard

E
Charging base
Ladestation
Station de recharge
Base di ricarica
Oplaadstation

F
Charging travel case
Box für Reiseladegerät
Coffret - chargeur de voyage
Custodia da viaggio
Oplaadreisetui

**BEFORE FIRST USE / VOR DEM ERSTEN GEBRAUCH /
AVANT LA PREMIÈRE UTILISATION / PRIMO UTILIZZO /
VOOR HET EERSTE GEBRUIK**



Refer to your user manual for additional instructions.





Nur bestimmte Gerätetypen Vor dem ersten Gebrauch Weitere Informationen finden Sie in der Bedienungsanleitung.

Reportez-vous au mode d'emploi pour de plus amples instructions.







Per ulteriori informazioni, consultate il manuale dell'utente.

Raadpleeg uw gebruiksaanwijzing voor extra instructies.


**USING / DAS GERÄT VERWENDEN / UTILISATION /
UTILIZZO / GEBRUIK**



**TRAVELING WITH YOUR SONICARE / REISEN MIT IHRER
SONICARE / VOYAGER AVEC VOTRE SONICARE /
IN VIAGGIO CON IL VOSTRO SONICARE / REIZEN
MET UW SONICARE**



**REPLACEMENT BRUSH HEADS / ERSATZBÜRSTENKÖPFE /
TÊTES DE BROSE DE RECHANGE / TESTINE DI RICAMBIO
PER SPAZZOLINO / VERVANGENDE OPZETBORSTELS**



Replace brush head every 3 months for optimal results.

Um optimale Ergebnisse zu erzielen, tauschen Sie den Bürstenkopf nach 3 Monaten aus.

Pour un résultat optimal, remplacez la tête de brosse tous les 3 mois.

Per ottenere ottimi risultati, sostituite la testina ogni 3 mesi.

Vervang de opzetborstel om de 3 maanden voor optimale resultaten.



4235.020.4914.1

www.philips.com/Sonicare