



Philips  
Combi air purifier and  
humidifier

**CADR 246 m<sup>3</sup>/h**

Formaldehyde removal 90% (3h)  
PM2.5 removal rate >99%

**AC4081**

# Healthy air always

With smart air and humidity control

Breathe healthy air with the Philips combi air purifier and humidifier. Its unique VitaShield IPS technology cleans while the humidifier moistens the air comfortably. The Healthy air protect alert and lock assure you of healthy air always.

## **Always**

- Child lock prevents unintended changes in settings
- Healthy air protect alert/lock

## **Healthy Air**

- Smart air and humidity control purifies and moistens air
- Dual sensors give clear feedback on air quality and humidity
- Superior VitaShield IPS with German technology

## **Simplicity**

- 3-step humidity setting
- 5 settings to give you ample choice in fan speed.
- 3-step LED indicators show air quality level clearly
- Auto mode ensures hassle-free functioning.
- lower speed & noise level in sleep mode with all lights off
- 1/4/8 hour easy-to-set timer

# PHILIPS

# Specifications

## Design specifications

- Control panel: Button
- Fan speed indicators: On control panel
- Materials of main body: ABS plastic
- Color(s): Chocolate brown

## Technical specifications

- Voltage: 110 V
- Frequency: 60 Hz
- Cord length: 1.8 m
- Air quality sensor(s): Particle
- Wattage (China): 47.5 W
- Noise level: 31.8 - 57.2 dB(A)

## Weight and dimensions

- Product weight: 8.3 kg
- F-box weight (incl. product): 10.2 kg

- Product dimensions (W x D x H): 210 x 365 x 605 mm
- F-box dimensions (W x D x H): 291 x 447 x 672 mm

## Logistic data

- 12 NC (China): 8834 081 80660
- Country of origin: China

## Performance

- Filters out TVOCs: 99%\*\*\*
- CADR (cigarette smoke): 246 m³/h\*\*\*
- Filters out PM2.5: 99%\*\*\* in 40 mins
- Filters out toluene: 99%\*\*\* in 1hr
- Filters out viruses: 99.1%\*\*\*\*
- Filters out bacteria: 99%\*\*\*\*
- Humidification rate: 339 ml/hr

# Highlights

## Dual sensors

Dual sensors measure and control air quality and humidity and give clear feedback on air quality and humidity

## Healthy air

Healthy air protect alert gives you a timely warning when to replace the filter. When the filter is not replaced in time, the appliance stops functioning to avoid running with no effect wince filter is full. It is your assurance of healthy air always.

## Smart air and humidity control

Smart air and humidity control purifies and moistens air

## Child lock

Child lock prevents unintended changes in settings.

## Air quality indicator

3-step LED indicators show air quality level clearly

## Auto mode

Auto mode ensures hassle-free functioning.

## Sleep mode

Product turns down fan speed and noise level in sleep mode, all lights on products are off.

## 3-step humidity setting

3-step humidity setting adjusts humidity to your liking


## 5 settings in fan speed

5 settings to give you ample choice in fan speed.

## Timer

1/4/8 hour easy-to-set timer

## Superior VitaShield IPS



VitaShield IPS

Superior VitaShield IPS with German technology effectively filters out harmful gases (formaldehyde, toluene and TVOC etc.) and pollutants ( ultrafine particles, allergens, bacteria and viruses) with a size of >0.02 micron. \*



Issue date 2016-02-14

Version: 1.0.2

© 2016 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* \*\*CADR (Cigarette) is tested by third party \*\*\*PM2.5 is tested by third party (5.0 +/- 0.5)mg \*\*\*Formaldehyde, Toluene and TVOC are tested by third party in 30 M3 lab in 3 hours (1.0 +/- 0.2)mg/m³, (2.0 +/- 0.4)mg/m³ and (6.0 +/- 1.2)mg/m³ (all according to GB/T 18801-2008, GB/T18883-2002)

\* \*According to 2008 Microbiological Risk Assessment Report of the World Health Organization (WHO), the avian influenza, human influenza viruses, Legionella, Hepatitis viruses and SARS coronavirus are larger than 0.02 micron.

\* \*\*\*Anti-bacterial properties tested according to GB21551.3 using Staphylococcus Albus, initial concentration is 1\*10<sup>5</sup> cfu/m³

\*\*\*\*Tested for MS2 virus species removal at Intertek Institute, USA (tested 2hrs in 11 m³ chamber), initial concentration is 1.3\*10<sup>7</sup> pfu/ml)

\* Tested for Staphylococcus Albus and MS2 virus, condition refer to \*\*\*\*